








SP-12&20

No. Catalog	Product Name	Specifications	Appearance	Description
SP-12	Liquid-Based Cytology Slide Preparation System	12 slides		The system can remove the impurities that will affect diagnosis by centrifugation and sedimentation to get cells with diagnostic value.
SP-20		20 slides		

Consumables for SP-12&20

No. Catalog	Product Name	Specifications	Appearance	Description
LC-01-01	Cell Preservation Solution	15ml/bottle		The cells can be preserved in good condition with the solution.
LC-02-01	Cytology slide preparation device	14mm		The target cells can be fixed very well on the slide.
LC-04-01	Sedimentation buffer solution	500ml/bottle		Make the cells separately by gravity.

Consumables for SP-M2

No. Catalog	Product Name	Specifications	Appearance	Description
LC-01-01	Cell Preservation Solution	15ml/bottle		The cells can be preserved in good condition with the solution.
LC-05-01	Membrane Filter	14mm		The target cells can be fixed very well on the slide.
LC-05-01	Adhesive Microscope Slides	50pcs/kit		Adhesive Microscope Slides

SEEKARY

Hangzhou Testsea Biotechnology Co.,Ltd
No. 13-2 Guanshan Road, Yuhang District, Hangzhou, 311115, China

Kevin.dong@testsealabs.com
+86 13003611301
www.testsealabs.com



TESTSEALABS®